




**E-Charm
Ceramic
Self-Ligating
Brackets**

E-Charm Self-ligating Brackets



Reliable bonding base design enhance bonding force


Manufactured by ceramic injection molding provides smooth, rounded surface and corners design

No metal parts showing on the door

Low profile design, improve the comfort of patients

Tie wing area satisfies different traction and ligation

■ Identification System



E-Charm Self-ligating Bracket

Item No. : EC-201/EC-202

	Bracket	Torq.	Ang.	Width(mm)	0.022R	0.022L	System
Maxillary	Central	+12°	+3°	3.20	EC-201-211	EC-201-221	Roth
	Central	+17°	+4°	3.20	EC-202-211	EC-202-221	MBT
	Lateral	+ 5°	+8°	3.00	EC-201-212	EC-201-222	Roth
	Lateral	+10°	+8°	3.00	EC-202-212	EC-202-222	MBT
	Cuspid	0°	+8°	3.00	EC-201-213	EC-201-223	Roth/MBT
	1st Bicuspid	- 7°	0°	3.00	EC-201-214	EC-201-224	Roth/MBT
	2nd Bicuspid	- 7°	0°	3.00	EC-201-215	EC-201-225	Roth/MBT
Mandibular	Central	- 1°	0°	2.80	EC-201-241	EC-201-231	Roth
	Central	- 6°	0°	2.80	EC-202-241	EC-202-231	MBT
	Lateral	- 1°	0°	2.80	EC-201-242	EC-201-232	Roth
	Lateral	- 6°	0°	2.80	EC-201-242	EC-201-232	MBT
	Cuspid	- 6°	+3°	3.00	EC-201-243	EC-201-233	Roth/MBT
	-	-	-	-	-	-	-
	-	-	-	-	-	-	-

System	Slot	Full Set Item Number
Roth	0.022	EC-201-2B (3# with hook)
MBT	0.022	EC-202-2B (3# with hook)

