

ReliaBox Quick Start Guide



# Unlock the benefits of ReliaBox

The key to efficient pre-pasted bracket placement



**ReliaBox innovative storage solution is designed to revolutionize your bracket application process.**

Pre-paste brackets for patient appointments up to a week in advance, reducing chairside time and increasing convenience for you and your team.



### ✓ Smart & Systematic Design.

ReliaBox is divided into 4 quadrants with 4 anterior tweezers and 3 posterior tweezers per quadrant.

### ✓ Accessible & Versatile.

When hung correctly, tweezers protrude 1-2" for convenient chairside access and hangers are designed to be used with any bracket holders, not just the Reliance brand.

### ✓ Protects Against Ambient Air & Light.

The custom designed locking seal allows for clinicians to pre-paste brackets up to 7 days in advance - keeping out both ambient air and light.

### ✓ Efficient & Effective Chairside.

Shave at least 15 minutes off every chairside appointment utilizing the ReliaBox pre-paste system!  
\*IMPORTANT: Do not place used tweezers back in box. Place with other used instruments for sterilization.



See it in action



NEW PRODUCT TIPS TO:

# Unlock the benefits of ReliaBox

[ReliaBox Video  
Link](#)



# ReliaBox

The key to efficient pre-pasted bracket placement



31.5cm

20.5cm

13cm



# ReliaBox

The key to efficient pre-pasted bracket placement

Ref. RBX: **ReliaBox** with spaces for 28 reverse-action tweezers for pre-coating/pre-pasting Brackets. No Tweezers included in this item number.



Ref. RBXK: **ReliaBox** with 16 anterior reverse-action tweezers and 12 posterior reverse-action tweezers.  
(Not for sale to EU, CH, UK, nor Australia)



# ReliaBox

The key to efficient pre-pasted bracket placement



RBTA Anterior

RBTA 1 x Anterior reverse action tweezers for use with ReliaBox. (Not for sale to EU, CH, UK, nor Australia)



RBTP Posterior

RBTP 1 x posterior reverse action tweezers for use with ReliaBox. (Not for sale to EU, CH, UK, nor Australia)

RBTA and RBTP Sold Individually

# ReliaBox

The key to efficient pre-pasted bracket placement

- **Freedom to pre-coat / pre-paste any bracket brand.**

Ortho professionals can now pre-coat their own brackets, and buccal tubes whatever the brand with GoTo, They are no longer limited to one pre-coated bracket or buccal-tube brand,

- **Pre-pasting up to 7 days before bonding**

The pre-coated / pre-pasted case of brackets and buccal tubes in the **ReliaBox** can be stored ready-to-use for up to 7 days.

- **Now available for international sale**



**RELIANCE**®

# ReliaBox

The key to efficient pre-pasted bracket placement

## - Cost-efficient

The cost of the box itself pays for itself in the two patient bondings with the savings made in chairside time and far-reduced costs compared with pre-coated brackets currently available.

## - Time-saving

Approx. 45% less time spent on bonding. Approx. 40% less chair-time, leaving orthodontists more time to:

1. spend with each patient and
2. to see more patients.



RELIANCE®



# ReliaBox

The key to efficient pre-pasted bracket placement



# GoTo<sup>®</sup>

Light Cure Bracket Adhesive

ReliaBox and GoTo – The perfect team !!!



 **RELIANCE**<sup>®</sup>

RELIANCE